

Auria® Harmony™ Sound Processor

Compatibility Chart



- ✓ Auria T-Mic® Microphone
- ✓ Auria Standard Earhook
- ✓ Auria Direct Connect Earhook
- ✓ Auria Direct Connect Cable
- ✓ Audio Interface Cable
- ✓ Auria T-Coil
- ✓ Auria iConnect
- ✓ Auria PowerCel™ Slim
- ✓ Auria PowerCel Plus
- ✓ Auria PowerCel Charging System
- ✓ Battery Charger Car Adaptor
- ✓ Auria Sport Case
- ✓ Auria Carrying Case
- ✓ Auria PowerCel Key Chain Wallet
- ✓ Auria Accent Color Covers and Headpiece Caps
- ✓ Little Things Wallet

Please note that the Harmony Processor (CI-5235) is NOT COMPATIBLE with:



Auria Firefly® (CI-5720)



Auria AA PowerPak™ (CI-7400)



System Sensor (CI-5820)

- Using the Firefly and System Sensor with the Harmony Processor may result in false diagnostics.
- If the AA PowerPak is used to power the Harmony Processor, the processor's LED Status Indicator may incorrectly indicate a low battery.
- In rare instances, using the AA PowerPak with the Harmony Processor may result in tripping the PowerPak's fuse.

If you have any questions, please contact Customer Care at 1(877) 829-0026 or customerservice@advancedbionics.com.